


AFKa Bench and Floor Weighing Scale



Speed, Performance, Value - Perfect Balance

- 6 Weighing units (kg, g, lb, oz, lb:oz, Newtons)
- Colour coded keys
- Hold function
- Large LCD display
- Stainless steel weighing pan
- Check weighing LED bar graph
- External calibration
- RS-232 Output
- Rechargeable battery up to 90 hours operation
- Weighing
- Parts counting
- Check weighing
- Animal / dynamic weighing
- Percentage weighing



Technical Specifications


Model	AFK 165a	AFK 330a	AFK 660a	AFK 1320a
Maximum capacity	165lb/75kg	330lb/150kg	660lb/300kg	1320lb/600kg
Readability	0.01lb/5g	0.02lb/10g	0.05lb/20g	0.1lb/50g
Repeatability (S.D.)	5g	10g	20g	50g
Pan size	15.7 x 19.7" / 400 x 500mm			23.6 x 31.5" / 600 x 800mm

General Specifications

Typical stabilisation time	2-3 Sec
Units of measurement	Gram (g), Kilogram (kg), Pound (lb), Pound-Ounce (lb-oz), Newton (N)
Operating temperature	32°F-104°F
Power supply	12vDC 800mA adapter or rechargeable battery (6v 4.5Ah battery)
Display	Backlit LCD Display with 1.6" / 40mm digits and capacity tracker
Dimensions notes	AFK165a-660a 15.7" x 24.4" x 31.1", AFK1320a 23.6

Key Accessories

- In-use Cover
- RS-232 Cable to PC
- RS-232 to USB Interface Cable (must be ordered in conjunction with the scale)
- ATP Thermal Printer
- Thermal Paper For Adam ATP (10 pack)
- Adam DU - Data Collection Program
- AE Relay Box



Adam Equipment - offices worldwide

www.adamequipment.com