

# Replacement Consumables for Millipore Lab Water Systems

Rephile offers a wide range of replacement cartridges and accessories for Millipore lab water systems.

## New offerings



with RFID tags

**RephiDuo P Pack** compatible to Millipore PROGOTOS2

**Designed for water softening and to remove large particles, organic compounds and chlorine from tap water. The Cat. No. is RR800CP01.**

Compatible with the following Millipore water systems:

- Elix Advantage 3/5/10/15
- Elix Reference 3/5/10/15
- Milli-Q Integral 3/5/10/15



with RFID tags

**RephiDuo U Pack** compatible to Millipore Q-Gard and QPAK

**Designed to remove ionic and organic contaminants down to trace levels from RO or distilled pure water.**

Compatible with the following Millipore water systems:


- Milli-Q Advantage A10
- Milli-Q Reference
- Milli-Q Direct 8/16

## Repair services for Millipore PCB and power boards

**Rephile provides repair services for failed control boards and power boards of various Millipore water systems and warrant full functionality of repaired boards for 6 months after service.**

Millipore Cat. No.	Description	Millipore systems
ZF3000022	PCB ELIX/RiOs INTERFACE	RiOs (96 / Century) Elix (96 / century)
ZF3000023	PCB MILLI-Q INTERFACE	Milli-Q (96 / Century)
ZF3000405 ZF3000487	MQ ADV_REF PWR BOARD	Milli-Q Advantage A10 Milli-Q Reference
ZF3000406	MQ ADV MAIN CONTROL BOARD	Milli-Q Integral (3/5/10/15) Milli-Q Direct (8/16) Milli-Q Advantage A10 Milli-Q Reference Elix Advantage (3/5/10/15) Elix Reference (3/5/10/15)
ZF3000433	MQ INT POWER SUPPLY BOARD	Milli-Q Integral (3/5/10/15) Milli-Q Direct (8/16) Elix Advantage (3/5/10/15) Elix Reference (3/5/10/15)

For any further information, please contact us by [info@rephile.com](mailto:info@rephile.com).



**A10 UV Lamp**



**Vent Filter**



**EDI**



**Final Filter**

## Cartridges and Filters

Millipore System	Millipore Cat. No.	RephiLe Compatible Cat. No.	Product Description
Milli-Q Academic Milli-Q Biocel Milli-Q Gradient Milli-Q Synthesis Milli-Q Element	QGARD00R1	RR400Q101	RephiDuo U Pack 1, compatible to QGARD00R1
	QGARD00D2	RR400Q201	RephiDuo U Pack 2, compatible to QGARD00D2
	QTUM000IX	RRQTUM0IX	RephiSolo U Pack IX, compatible to QTUM000IX
	QTUM000EX	RRQTUM0EX	RephiSolo U Pack EX, compatible to QTUM000EX
	QTUMMPEIX	RRQTUMFIX	RephiSolo U Pack IX kit, compatible to QTUMMPEIX
Milli-Q Advantage A10 Milli-Q Reference	QGARDT1X1	RR800Q101	RephiDuo U Pack T1, compatible to QGARDT1X1
	QGARDT2X1	RR800Q201	RephiDuo U Pack T2, compatible to QGARDT2X1
Milli-Q Advantage A10 Milli-Q Integral Milli-Q Reference	QTUM0TIX1	M1001MTIX	RephiSolo U Pack TIX, compatible to QTUM0TIX1
	QTUM0TEX1	M1001MTEX	RephiSolo U Pack TEX, compatible to QTUM0TEX1
Milli-Q Direct	PROG000T3	M1001CPT3	RephiSolo P Pack T3, compatible to PROG000T3
	QPAK00TIX	RR800QTIX	RephiDuo U Pack TIX, compatible to QPAK00TIX
	QPAK00TEX	RR800QTEX	RephiDuo U Pack TEX, compatible to QPAK00TEX
Milli-Q Plus Milli-Q PF Plus Milli-Q UF Plus Milli-Q UV Plus	CPMQ004R1	RR100Q101	RephiQuatro U Pack 1, compatible to CPMQ004R1
	CPMQ005R1	RR100Q1KT	RephiQuatro U Pack 1 kit, compatible to CPMQ005R1
	CPMQ004D2	RR100Q201	RephiQuatro U Pack 2, compatible to CPMQ004D2
	CPMQ005D2	RR100Q2KT	RephiQuatro U Pack 2 kit, compatible to CPMQ005D2
Elix Advantage Elix Reference Milli-Q Integral	PROG0T0S2	RR800CP01	RephiDuo P Pack S2, compatible to PROG0T0S2
Elix 3/5/10 RiOs 5/8/16 Direct-Q 5	PROG00002	RR400CP01	RephiDuo P Pack 2, compatible to PROG00002
	PROG0NP02	RR400NP01	RephiDuo P Pack 2 for softer water, compatible to PROG0NP02
	PROG000S2	RR400CPS2	RephiDuo P Pack S2 (bacteriostatic activated carbon), compatible to PROG000S2
	PROG00001	RR400CP00	RephiDuo P Pack 1, compatible to PROG00001, 2/pk
Elix Clinical Elix Gulfstream Elix 20/35/70/100 RiOs 30/50/100/150/200	PROGTL0S1	RR500CP01	RephiDuo P Pack TL1, compatible to PROGTL0S1
	PROGTLCS1	RR500CPC1	RephiDuo P Pack TL1 with RO cleaning media, compatible to PROGTL0S1
	PRPK0L0S1	RR500CPP1	RephiDuo P Pack L1, compatible to PRPK0L0S1
	QGARDTL01	RR500Q101	RephiDuo U Pack TL, compatible to QGARDTL01
Simplicity Simplicity 185 Synergy Synergy 185	SIMPAK0R1	RRSYNM0IX	RephiDuo S Pack 1, compatible to SIMPAK0R1
	SIMPAK0R1	RRSYNMFIX	RephiDuo S Pack 1 kit, compatible to SIMPAK0R1
	SIMPAK0D2	RRSYNM0DX	RephiDuo S Pack 2, compatible to SIMPAK0D2
	SIMPAK0D2	RRSYNMF0DX	RephiDuo S Pack 2 kit, compatible to SIMPAK0D2
Simplicity(SMART) Simplicity UV	SIPK0SIA1	M1001SIA1	RephiSolo SI Pack 1, compatible to SIPK0SIA1
	SIPK0SIX2	M1001SIX2	RephiSolo SI Pack 2, compatible to SIPK0SIX2
Synergy(SMART) Synergy UV	SYPK0SIA1	M1002SIA1	RephiDuo SI Pack 1, compatible to SYPK0SIA1
	SYPK0SIX2	M1002SIX2	RephiDuo SI Pack 2, compatible to SYPK0SIX2
AFS 10/10D/18/18D/30/30D AFS-25 EDI/AFS-50 EDI Milli-RO 3/5/8/12 Milli-RO 30/60/90 Milli-RX 20/45/75	CPR0P0402	RR100CP01	RephiQuatro P Pack, compatible to CPR0P0402
	CPR0NP402	RR100NP01	RephiQuatro P Pack for softer water, compatible to CPR0NP402
AFS 3/8/16	PROG00003	RR400CP03	RephiSolo P Pack 3, compatible to PROG00003
	CP1ALLRES	RR1ALLRES	RephiSolo Clinical Pack, compatible to CP1ALLRES
	CP2ALLRES	RR400Q401	RephiDuo Clinical Pack, compatible to CP2ALLRES
AFS 10/18/30 AFS 10D/18D/30D AFS-25 EDI/AFS-50 EDI	CP4ALLRES	RR100Q401	RephiQuatro Clinical Pack, compatible to CP4ALLRES
Milli-DI	CPDI000L1	RRDI000L1	RephiDuo DI Pack L, compatible to CPDI000L1
	PRPK00001	RRPRPK001	RephiDuo P Pack pretreatment, compatible to PRPK00001, 2 per/box
	ZLFI00001	RASDI0001	RephiLe SDI Kit
Elix Milli-Q Integral	TANKMPK01	RATANKVN1	Tank Vent Filter with CO <sub>2</sub> scavenger, compatible to TANKMPK01
RiOs Milli-Q Direct	TANKMPK02	RATANKVN2	Tank Vent Filter, compatible to TANKMPK02
Simplicity	SLFH02510	RJF324510	Vent Filter, compatible to SLFH02510, 10/pk
Ultrapure Water Systems POD	MPGP04001 MPGP02001	RAFFC0250	0.2 µm PES High Flus Capsule Filter, compatible to MPGP04001, MPGP02001
Simplicity	MPGL04001 SIMFILTER	RAFFC0200	0.2 µm PES Capsule Filter, compatible to MPGL04001, SIMFILTER

## UV Lamps, RO Membranes, Pumps and Others

Millipore System	Millipore Cat. No.	RephiLe Compatible Cat. No.	Product Description
AFS Direct-Q 3/5 Elix (UV) / RiOs Elix Advantage / Elix Reference Elix Essential / RiOs Essential Milli-Q Direct Milli-Q Integral Milli-RO 3/5	CDRC00201	RAR00201	Small RO membrane, compatible to CDRC00201
	CDRC60201	RAR060201	RO membrane for small Elix and RiOs, compatible to CDRC60201
	CDRC60202	RAR060202	RO membranes for small Elix and RiOs, compatible to CDRC60202, 2/pk
	CDRC75202	RAR075202	RO membranes for Elix 15, compatible to CDRC75202, 2/pk
Elix 35/70/100 RiOs 50/100/150/200	CDRC050HG	RAR0050HG (Reuse existing housing)	RO membrane w/o housing, compatible to CDRC050HG
Milli-RO 6/8/10/12	CDRC012HG	RAR0012HG (Reuse existing housing)	RO membrane w/o housing, compatible to CDRC012HG
Milli-RX 20/45/70 Milli-RO 30/60/90	CDRCH12HG	RAR0H12HG (Reuse existing housing)	RO membrane w/o housing, compatible to CDRCH12HG
Milli-Q Gradient (Century) Milli-Q Synthesis (Century) Milli-Q Element Milli-Q Advantage A10 Milli-Q Direct Milli-Q Integral Milli-Q Reference	ZMQUVP01	265021701	UV Lamp 185/254 nm, compatible to ZMQUVP01
Milli-Q Gradient (96) Milli-Q Synthesis (96)	QUVLQSL01	RAUVLQSL1	Replacement UV Lamp, compatible to QUVLQSL01
Direct-Q 3 UV Simplicity UV (SMART) Synergy UV (SMART)	SYN185UV1	265020601	UV Lamp 185 nm, compatible to SYN185UV1
Synergy UV (High Flow)	ZLUVWPM01	RAUV212B1	UV Lamp 185 nm, compatible to ZLUVWPM01
Milli-Q (96/Century) A10 Milli-Q Advantage A10 Milli-Q Integral	ZFA10UVM1 ZFA10UV01	RAA10UVM1	UV Lamp for A10 TOC monitor, compatible to ZFA10UVM1, ZFA10UV01
Elix UV Elix Advantage/Elix Reference Elix Essential Milli-Q Direct Milli-Q Integral RiOs-DI 3 UV	ZLXUVP01	265010601	UV Lamp 254 nm, compatible to ZLXUVP01
Elix 20/35/70/100 RiOs 30/50/100/150/200 ASM (Cat. No.: TANKASMIN)	ZLXUVP01	265011001	UV Lamp 254 nm, compatible to ZLXUVP01
ASM (Cat. No.: TANKS50UV/ TANKS5LUV/ZFRES50UV/ZFRES60UV)	ZFRES00UV	2650110C1	UV Lamp 254 nm, compatible to ZFRES00UV
Elix (UV) Elix Advantage Elix Reference Elix Essential Milli-Q Integral AFS 10E/15E	ZLX0EDI03	W3T101571	LPX 03 CEDI Replacement Module, compatible to ZLX0EDI03
	ZLX0EDI05	W3T101572	LPX 05 CEDI Replacement Module, compatible to ZLX0EDI05
	ZLX0EDI10	W3T101573	LPX 10 CEDI Replacement Module, compatible to ZLX0EDI10
	ZLX0EDI15	W3T262701	LPX 15 CEDI Replacement Module, compatible to ZLX0EDI15
30 PE reservoirs	FTPF05380	RAPF05380	30 L Tank Level Sensor, compatible to FTFP05380
60 PE reservoirs	FTPF05381	RAPF05381	60 L Tank Level Sensor, compatible to FTFP05381
100 PE reservoirs	FTPF06805	RAPF06805	100 L Tank Level Sensor, compatible to FTFP06805
SDS 200	FTPF07743	RAPF07743	200 L Tank Level Sensor, compatible to FTFP07743
	ZF2000022	RAPF20022	200 L Tank Level Sensor, compatible to ZF2000022
SDS 350	FTPF07744	RAPF07744	350 L Tank Level Sensor, compatible to FTFP07744
	ZF2000023	RAPF20023	350 L Tank Level Sensor, compatible to ZF2000023
Elix 3/5/10 RiOs 5/8/16	ZF3000000	RASP00000 (Reuse existing housing)	Small RiOs / Elix pump, 36VDC w/o support, compatible to ZF3000000
Milli-Q (96/Century)	ZF3000001	RASP00001 (Reuse existing housing)	Milli-Q 96 / Century pump, 36VDC w/o support, compatible to ZF3000001
Elix Advantage Elix Reference Milli-Q Integral Milli-Q Direct	ZF3000431	RASP00431 (Reuse existing housing)	RO Booster Pump compatible to Millipore ZF3000431 (RephiLe Cat. RASP00431)
Milli-Q Advantage Milli-Q Direct Milli-Q Integral Milli-Q Reference	ZF3000401	RASP00401 (Reuse existing housing)	Recirculation Pump compatible to Millipore ZF3000401 (RephiLe Cat. RASP00401)
Milli-Q Synergy	ZFMQ000PR	RAPR58561	Pressure Regulator, compatible to ZFMQ000PR

## Ultrapure Cartridges

- Filled with LeFil™ and OrganeFil™ media, which are composed of ion exchange resins and activated carbon
- Designed to remove specific types of water contaminants with large capacity



### RephiSolo and RephiDuo U Pack

Compatible to Millipore Q-Gard & Quantum

- Compatible with the following Millipore water systems:

- Milli-Q Academic
- Milli-Q Bioceel
- Milli-Q Gradient
- Milli-Q Synthesis
- Milli-Q Element



### RephiQuatro U Pack

Compatible to Millipore QPAK

- Compatible with the following Millipore water systems:

- Milli-Q Plus
- Milli-Q PF Plus
- Milli-Q UF Plus
- Milli-Q UV Plus
- AFS 200
- AFS 300
- Alpha-Q
- Milli-Q RG



### RephiDuo S Pack

Compatible to Millipore SimPak

- Compatible with the following Millipore water systems:

- Simplicity
- Synergy

## Pretreatment Filtration Cartridges

- Designed to remove most particles, organic compounds and chlorine from tap water to protect the RO membranes

### RephiDuo P Pack

Contains activated carbon and depth filtration media in a polypropylene housing.



- Compatible to Millipore Progard
- Compatible with the following Millipore water systems:
  - Elix 3/5/10
  - RiOs 5/8/16
  - Direct-Q 5

### RephiQuatro P Pack

Designed in a duplex configuration, these cartridges can be rotated and reconnected when one-half has been exhausted.



- Compatible to Millipore ROPAK
- Compatible with the following Millipore water systems:
  - Milli-RX 20/45/75
  - Milli-RO 5/8/12/30/60/90
  - Milli-U 6/10

### RephiDuo P Pack TL1

Suitable for large water systems and designed to ensure long service time.



- Compatible to Millipore Progard TL
- Compatible with the following Millipore water systems:
  - Elix 20/35/70/100
  - RiOs 30/50/100/150/200
  - Elix Clinical
  - Elix Gulfstream Clinical

Superior  
Quality

Guaranteed  
Performance

Highly  
Cost Effective

All rights reserved © Copyright 2016 Rephile Bioscience, Ltd.

**Rephile Bioscience, Ltd.**

Toll Free: +1-855-rephile (+1-855-737-4453)

Email: [info@rephile.com](mailto:info@rephile.com)

Lit. No.: RFPR0871409