

Speed, Performance, Value - Perfect Balance

- 3 Weighing units
- Full range tare
- Backlit LCD display
- Stainless steel pan
- ShockProtect™ overload three-point protection
- Lock down mounting slot for Kensington® type lock and cable
- Removable draught shield (120mm pan size only)
- Lightweight and portable
- Modern low profile design
- Colour coded, sealed keypad
- External calibration
- Auto sleep / power down function to save battery life
- AC Adapter
- Low battery indication
- Below balance weighing with hanger
- Weighing
- Below balance weighing



Technical Specifications

Maximum capacity	1750g
Readability	0.05g
Repeatability (S.D.)	0.1g
Linearity (\pm)	0.1g
Pan Size	4.7" Ø / 120mm Ø

General Specifications

Model	CQT 1752GR
Typical stabilisation time	2-3 Sec
Power supply	12VAC 150mA

Key Accessories

- M1 1kg Calibration Weight
- Hard Carry Case with Lock
- Security Lock and Cable
- In-use Cover
- ASTM 2 - 1000g calibration weight



Adam Equipment - offices worldwide

www.adamequipment.com