




# MDW-160M Mechanical Person Scale



## Speed, Performance, Value - Perfect Balance

- Height rod as standard (measures up to 60-212 cm)
- Robust metal housing
- Anti-slip platform
- Adjustable zero point
- Wheels to allow easy movement
- 2 weighing units (kg, lb)
- Weighing
- Measuring height



## Technical Specifications

Maximum capacity	350lb/160kg
Readability	4oz/100g
Repeatability (S.D.)	0.4lb / 200g
Linearity (±)	0.4lb / 200g
Platform size	14.8 x 10.8" / 375 x 275mm

## General Specifications

Model	MDW-160M
Units of measurement	Kilogram (kg), Pound (lb)
Division of measurement	0.1cm
Housing	Metal
Overall dimensions (wxdxh)	20.9 x 10.8 x 58.5" / 530 x 275 x 1485mm



Adam Equipment - offices worldwide

[www.adamequipment.com](http://www.adamequipment.com)