



Speed, Performance, Value - Perfect Balance

- Ultra sonic sensor measures height up to 210cm and automatically calculates BMI
- Nonslip Platform
- Backlit LCD display
- BMI-Body mass index
- Rechargeable battery
- Weighing units include lbs, lbs/ozs and kg
- Security locking point to secure the scale in one place
- Measures height in feet, inches and centimetres
- Weighing
- Measuring height



Technical Specifications

Maximum capacity	660lb/300kg
Readability	0.1lb/50g
Repeatability (S.D.)	50g

General Specifications

Model	MUW 300L
Typical stabilisation time	2-3 Sec
Power supply	12VDC @500mA AC adapter, internal rechargeable battery
Display	Backlit LCD with 3 lines for weight, height and BMI
Overall dimensions (wxdxh)	91.3 x 14.8 x 23.6" / 2320 x 375 x 600mm

Key Accessories

- RS-232 to USB Interface Cable (must be ordered in conjunction with the a
- ATP Thermal Printer
- Thermal Paper For Adam ATP (10 pack)
- RS-232 to PC Cable



S232 cable)



Adam Equipment - offices worldwide

www.adamequipment.com