



Speed, Performance, Value - Perfect Balance

- 3 Weighing units (g, lb, oz)
- Hold function
- Full range tare
- Simple user-friendly operation
- Large easy to read LCD display
- Simple 2 button operation
- Solid 1 piece construction is lightweight and easy to clean
- Easy grip carrying handles
- Splashproof to protect from accidental spills
- Lightweight and portable
- Low cost and simple to use
- External calibration
- Battery power (2x AA Batteries)
- Auto sleep / power down function to save battery life
- DC Adapter
- Weighing
- Animal / dynamic weighing



Technical Specifications

| | |
|----------------------|------------|
| Maximum capacity | 44lb/20kg |
| Readability | 0.01lb/10g |
| Repeatability (S.D.) | 10 |
| Linearity (\pm) | 20g |

General Specifications

| | |
|----------------------------|--|
| Model | MXB Baby |
| Units of measurement | Kilogram (kg), Pound (lb), Pound-Ounce (lb-oz) |
| Power supply | 6V AC adaptor or 4 x AA Batteries |
| Display | LCD |
| Overall dimensions (wxdxh) | 23.6 x 13.4 x 2.8" / 600 x 340 x 70mm |
| Net weight | 4.4lb |



Adam Equipment - offices worldwide

www.adamequipment.com