



Speed, Performance, Value - Perfect Balance

- Low cost and simple to use
- Magnetic damping for faster stability
- Rigid construction for durability
- Zero adjust facility to enable zero reading
- Notched and tiered graduated beams with centre reading
- Stainless steel weighing pan
- Ruler for measuring cm
- Weight holders for storing additional weights
- Tripod legs for density: TBB 2610S and TBB 2610T models only.
- Large stainless steel pan
- Lock down mounting slot for Kensington® 482; type lock and cable
- Weighing
- Below balance weighing
- Density of solids



Technical Specifications

Model	TBB 2610T	TBB 2610S	TBB 610S
Maximum capacity	2610g	2610g	610g
Readability	0.1g	0.1g	0.1g
Tare range	+225g	-	-
Repeatability (S.D.)	0.2g	0.2g	0.2g
Linearity (\pm)	0.2g	0.2g	0.2g
Pan size	6" Ø / 152mm Ø		
Net weight	5.5lb		3.8lb

General Specifications

Units of measurement	Gram (g)
Calibration	External
Display	Graduated Beam
Housing	Aluminium
Overall dimensions (wxdxh)	19.9 x 4.3 x 6.3" / 505 x 110 x 160mm

Key Accessories

- Security Lock and Cable
- Hard Carry Case with Lock



Adam Equipment - offices worldwide

www.adamequipment.com